NADAL[®] COVID-19 Ag+Influenza A/B plus Test (test cassette)

REF 243204N-20

de Gebrauchsanweisung	2	pt Instruções de Utilização	38
en Instructions for use	8	cs Návod k použití	44
fr Instructions d'utilisation	14	fi Käyttöohje	50
es Instrucciones de uso	20	sv Användarinstruktioner	56
it Istruzioni per l'uso	26	Symbols	63
pl Sposób użycia	32	Our Teams	64

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Carl-Zeiss-Strasse 12
47445 Moers
Germany

Moers
Tel: +49 (2841) 99820-0
Fax: +49 (2841) 99820-1

Regensburg
Tel: +49 941 29010-0
Fax: +49 941 29010-50

www.nal-vonminden.com
info@nal-vonminden.com

Directors:
Sandra von Minden
Roland Meißner
Thomas Zander

Commercial reg. Klevé
HRB 5679
Steuer-Nr. 244/133/00130
UST-ID-Nr. DE 189 016 086

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Rev. 2, 2021-01-18 OM

1. Intended Use

The NADAL® COVID-19 Ag+Influenza A/B plus Test is a lateral flow chromatographic immunoassay for the simultaneous qualitative detection of SARS-CoV-2 as well as influenza virus type A and B viral nucleoprotein antigens in human nasal, nasopharyngeal or oropharyngeal specimens (see section 12 'Limitations'). This test is intended for use as an aid in the diagnosis of infections with SARS-CoV-2 as well as influenza virus type A and B. Note that the concentration of viral nucleoprotein antigens may vary in the course of the disease and might fall below the detection limit of the test. Possible infectiousness of test subjects cannot be ruled out based on negative test results. The test procedure is not automated and requires no special training or qualification. The NADAL® COVID-19 Ag+Influenza A/B plus Test is designed for professional use only.

2. Introduction and Clinical Significance

COVID-19 (Corona Virus Disease) is the infectious disease caused by the recently discovered coronavirus SARS-CoV-2. The most common symptoms of COVID-19 are fever, dry cough, fatigue, sputum production, shortness of breath, sore throat and headache. Some patients may have myalgia, chills, nausea, nasal congestion and diarrhoea. These symptoms begin gradually and are mild in most of the cases. Some people become infected but do not develop any symptoms and do not feel unwell. Most people (about 80%) recover from the disease without special treatment. Approximately one in six people who get infected with COVID-19 becomes seriously ill and develops difficulty breathing. Elderly people, and those with pre-existing conditions, such as high blood pressure, heart problems or diabetes, are more likely to develop serious illness. COVID-19 is transmitted via respiratory droplets that are exhaled by infected people via coughing, sneezing or talking. Most estimates of the incubation period for COVID-19 range from 1 to 14 days, during which people might already be infectious without showing disease symptoms.

Influenza is a highly contagious viral infection of the upper respiratory tract which is characterised by antigen variability, seasonality and the impact on the general population. Of the two main types (A and B) of influenza viruses, influenza A subtypes are differentiated by the antigen variability of surface glycoproteins (haemagglutinin and neuraminidase). The influenza A virus is the most prevalent and is associated with the most serious epidemics. Influenza can cause severe complications such as bronchitis or pneumonia, particularly in children, the elderly or those with chronic respiratory diseases. However, it most commonly occurs as a mild viral infection transmitted by respiratory secretions through sneezing or coughing. There are many other viral infections that can mimic influenza symptoms, making laboratory tests necessary to distinguish it from other acute respiratory infections.

With a sensitivity and specificity of almost 100%, RT-PCR is the gold standard of laboratory diagnostics. However, as a laboratory method, it does not allow immediate diagnosis and requires qualified personnel and expensive laboratory equipment. Antiviral influenza drugs, which have been available since the late 1990s, are most effective if administered soon after the onset of the illness (within 48

hours). Rapid tests for the detection of influenza antigens can therefore enable successful antiviral therapy to be carried out in good time, thereby improving patient healthcare and reducing health costs.

3. Test Principle

The NADAL® COVID-19 Ag+Influenza A/B plus Test is a lateral flow chromatographic immunoassay for the qualitative detection of SARS-CoV-2 as well as influenza virus type A and B viral nucleoprotein antigens in human nasal, nasopharyngeal or oropharyngeal specimens.

The NADAL® COVID-19 Ag+Influenza A/B plus Test contains two internal test strips (left side of the test cassette: COVID-19 Ag Test; right side of the test cassette: Influenza A/B Test).

Anti-SARS-CoV-2 antibodies are immobilised in the test line region 'T' of the membrane. Anti-influenza virus type A and B antibodies are immobilised in the test line regions 'A' and 'B', respectively, of the membrane. A specimen is added to an extraction tube containing buffer in order to release SARS-CoV-2 as well as influenza virus type A and B viral nucleoprotein antigens. During the test, extracted antigens, if present in the specimen, bind to anti-SARS-CoV-2 antibodies and/or anti-influenza virus type A and/or B antibodies conjugated to coloured particles and precoated onto the sample pad of the test cassette. The mixture then migrates along the membrane chromatographically by capillary action and interacts with the reagents on the membrane. The complexes are then captured by anti-SARS-CoV-2 antibodies as well as by anti-influenza virus type A and B antibodies in the test line regions 'T' as well as 'A' and 'B', respectively. Excess coloured particles are captured in each control line region 'C'. The presence of a coloured line in the test line region 'T' and/or 'A' and/or 'B' indicates a positive result. The absence of a coloured line in the test line region 'T' and/or 'A' and/or 'B' indicates a negative result.

The formation of a coloured line in each control line region 'C' serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 20 NADAL® COVID-19 Ag+Influenza A/B plus test cassettes
- Additional material provided according to 93/42/EEC:
Due to possible supply shortages of COVID-19 related accessory medical products, the swab manufacturer might change. Therefore, the supplied swabs are from one of the manufacturers listed below.

a) 20 sterile swabs, CE 0197

Jiangsu Changfeng Medical Industry Co., Ltd
Touqiao Town, Guangling District, Yangzhou,
Jiangsu 225109 China (authorised EU
representative Llns Service & Consulting GmbH,
Obere Seegasse 34/2, 69124 Heidelberg,
Germany)

b) 20 sterile swabs, CE 0197

Jiangsu Rongye Technology Co., LTD, Touqiao
Town, Yangzhou City, Jiangsu Province, China
(authorised EU representative Riomavix S.L., Calle
de Almansa 55, 1D, Madrid 28039 Spain)

c) 20 sterile swabs, CE 0197

CITOTEST LABWARE MANUFACTURING CO., LTD
No.48, Xinxiu Road, Haimen, Jiangsu province
(authorised EU representative WellKang Ltd,
Enterprise Hub, NW Business Complex,
1 Beraghmore Rd., Derry, BT48 8SE, Northern
Ireland)

- 20 extraction tubes incl. dropper caps
- 20 buffer ampoules for single use (300 µL each)*
- 1 reagent holder
- 1 package insert

*Buffer containing the following preservative: ProClin™ 300:
<0.03%.

The detergents contained in the buffer lyse and neutralise the virus.

No hazard labelling is required according to Regulation (EC) № 1272/2008 CLP. Concentrations are below exemption threshold.

5. Additional Materials Required

- Timer

6. Storage & Stability

Test kits should be stored at 2-30°C until the indicated expiry date. Test cassettes are stable until the expiry date printed on the foil pouches. Test cassettes must remain in the sealed foil pouches until use. Do not freeze the test kit. Do not use tests beyond the expiry date indicated on the packaging. Care should be taken to protect test kit components from contamination. Do not use test kit components if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to inaccurate results.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the test procedure prior to testing.
- Do not use the test beyond the expiration date indicated on the packaging.
- Do not use test kit components if the primary packaging is damaged.
- Tests are for single use only.
- Do not add specimens to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new extraction tube for each specimen obtained.
- Do not substitute or mix components from different test kits.
- Do not use the buffer if it is discoloured or turbid. Discolouration or turbidity may be a sign of microbial contamination.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear appropriate personal protective equipment, such as face mask, isolation gown, gloves and eye protection during specimen collection, preparation and test procedure.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.

- Further specimen processing and patient management should follow local guidelines and regulations.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Temperature can adversely affect test results.
- Used testing materials should be disposed of according to local regulations.

8. Specimen Collection and Preparation

Nasal specimen:

- It is important to obtain as much secretion as possible. Insert the swab into the nostril.
- Gently push the swab until resistance is met at the level of the turbinate (up to 2.5 cm into the nostril).
- Rotate the swab 5 times against nasal wall to ensure that both mucus and cells are collected.
- Slowly withdraw the swab while continuing to rotate it.
- Repeat this process for the other nostril using the same swab in order to ensure that sufficient specimen volume is collected from both nasal cavities.

Nasopharyngeal specimen:

- Insert the swab into the nostril, parallel to the palate (not upwards) until resistance is encountered or the distance is equivalent to that from the ear to the nostril of the patient, indicating contact with the nasopharynx.
- Gently rub and roll the swab. Leave swab in place for several seconds to absorb secretions.
- Slowly remove the swab while rotating it. Specimens can be collected from both nostrils using the same swab, but it is not necessary to collect specimens from both sides if the tip is saturated with fluid from the first collection.

Oropharyngeal specimen:

- Gently insert a sterile swab into the pharynx and collect secretions by brushing the swab several times against the reddened posterior pharyngeal wall and both tonsillar pillars. Avoid touching the tongue, teeth and gums.

Note:

- Use only synthetic fibre swabs with plastic shafts. Do not use calcium alginate swabs or swabs with wooden shafts, as they may contain substances that inactivate some viruses and inhibit further testing.
- Swab specimens should be tested immediately after collection. Use freshly collected specimens for best test performance. If the test cannot be carried out immediately after specimen collection, the swab can be stored for up to a maximum of one hour in a clean, dry and sealed tube.
- Do not use specimens that are obviously contaminated with blood, as it may interfere with the flow of specimens and lead to inaccurate test results.

9. Test Procedure

Bring tests, specimens, buffer and/or controls to room temperature (15-30°C) prior to testing.

1. Place a clean extraction tube labeled with the patient or control identification into the designated area of the reagent holder.

2. Detach one buffer ampoule.

3. Open the ampoule by twisting the tip off.

4. Hold the buffer ampoule vertically over the tube to ensure that the entire buffer solution flows into the bottom part. Squeeze the buffer ampoule and add the entire buffer solution to the extraction tube without touching the edge of the tube.

5. Insert the swab with the collected specimen into the tube. Swirl the swab and squeeze it at least 5 times by compressing the wall of the extraction tube against the swab to extract the antigens contained in the swab.


6. Remove the swab, pressing it firmly against the wall of the tube to release as much liquid as possible. Dispose of the swab in accordance with guidelines for the handling of infectious agents.

Note: If the dropper cap is attached, the extracted specimen can be stored in the extraction tube for up to a maximum of 1 hour until the removal of the test cassette from its foil pouch.

7. Remove the test cassette from the foil pouch and use it as soon as possible. The best results will be obtained if the test is performed immediately after opening the foil pouch. Label the test cassette with the patient or control identification.

8. Place the test cassette on a clean and level surface.

9. Attach a dropper cap to the extraction tube, invert the tube and transfer 3 drops (approximately 60 µL) of the extracted solution to each specimen well (S) of the test cassette.

10. Start the timer.

11. Wait for the coloured line(s) to appear. Read the test result after 15 minutes. Do not interpret the result after more than 20 minutes.

15 min

10. Result Interpretation


Positive for COVID-19 Ag:

A coloured line develops in the control line region 'C' and another coloured line develops in the test line region 'T'.


Positive for influenza virus type A:

A coloured line develops in the control line region 'C' and another coloured line develops in the test line region 'A'.

Positive for influenza virus type B:

A coloured line develops in the control line region 'C' and another coloured line develops in the test line region 'B'.


Note: A combination of positive results for various parameters is possible.

Note: The colour intensity in the test line region 'T', 'A' or 'B' may vary depending on the concentration of SARS-CoV-2 as well as influenza virus type A and B viral nucleoprotein antigens in the specimen. Any shade of colour in the test line region 'T', 'A' or 'B' should be considered a positive result. Note that this is a qualitative test only and it cannot determine the analyte concentration in the specimen.


Negative:

A coloured line develops in each control line region 'C'. No lines develop in the test line regions 'T', 'A' and 'B'.

Invalid

The control line 'C' fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded. Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control line failure.

11. Quality Control

An internal procedural control is included in the test cassette:

A coloured line developing in each control line region 'C' is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Good laboratory practice (GLP) recommends the use of external control materials to ensure proper test kit performance.

12. Limitations

- The NADAL® COVID-19 Ag+Influenza A/B plus Test is for professional *in-vitro* diagnostic use only. The test should be used for the qualitative detection of SARS-CoV-2 as well as influenza virus type A and B viral nucleoprotein antigens in human nasal, nasopharyngeal or oropharyngeal specimens only. Neither the quantitative value nor the rate of increase/decrease in the concentration of SARS-CoV-2 as well as influenza virus type A and B viral nucleoprotein antigens can be determined with this qualitative test.
- The NADAL® COVID-19 Ag+Influenza A/B plus Test only detects the presence of SARS-CoV-2 as well as influenza virus type A and B viral nucleoprotein antigens in specimens and should not be used as the sole criterion for a diagnosis of COVID-19 and influenza A or B.
- Both viable and non-viable SARS-CoV-2 viruses as well as influenza viruses type A and B can be detected using the NADAL® COVID-19 Ag+Influenza A/B plus Test.
- The sections 'Specimen Collection and Preparation' as well as 'Test Procedure' must be followed closely while testing. Failure to follow them may lead to inaccurate test results because the antigen concentration in the swab is highly dependent on the correct procedure.
- As with all diagnostic tests, all results should be interpreted in conjunction with other clinical information available to the physician.
- In the course of SARS-CoV-2 or influenza A or B infection, the concentration of viral nucleoprotein antigens may fall below the detection limit of the test.
- If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of a SARS-CoV-2 or influenza A or B infection and should be confirmed via molecular assay.
- Additional testing is required to differentiate any specific influenza virus A and B subtypes or strains in consultation with state or local public health departments.
- Monoclonal antibodies may fail to detect or detect with less sensitivity influenza A virus antigens that have undergone minor amino acid changes in the target epitope region.
- Individuals who received nasally administered influenza A vaccine may have positive test results for up to three days after vaccination.
- Positive results do not preclude co-infections with other pathogens.
- Children tend to shed influenza virus more abundantly and for longer periods of time than adults. Therefore, testing specimens from adults may often yield lower sensitivity than testing specimens from children.
- Positive and negative predictive values are highly dependent on prevalence. The local prevalence should be taken into consideration when interpreting diagnostic test results.

13. Expected Values

SARS-CoV-2 viral particles are normally present in the respiratory tracts of COVID-19 patients. A positive test result can indicate an acute infection. Virus concentrations in nasal, nasopharyngeal or oropharyngeal swab specimens may vary in the course of the disease and might fall below the detection limit of rapid tests, even though patients are still showing

symptoms. Conversely, the virus might continue to be detectable over long periods of time even in convalescent patients. Possible infectiousness of test subjects cannot be ruled out based on negative test results.

On average, 3 to 10% of unvaccinated adults get infected with influenza each year (the number of infections varies significantly depending on the currently most prevalent strains). A positive rapid test result can indicate the presence of influenza virus type A and/or B viral nucleoprotein antigens in the specimen material, suggesting an acute infection with influenza A and/or B. Since patients might also acquire co-infections with other severe respiratory diseases, additional testing for other viral and bacterial infectious agents is recommended.

14. Performance Characteristics

Clinical performance

Diagnostic sensitivity and specificity

The NADAL® COVID-19 Ag+Influenza A/B plus Test was evaluated with clinical nasal swab specimens whose status was confirmed using RT-PCR.

The results are presented in the following tables.

		NADAL® COVID-19 Ag+Influenza A/B plus Test			
		Influenza A	Influenza B	SARS-CoV-2	
RT-PCR	No. of specimens				
Influenza A	Positive	42	38+/4-	0+/42-	0+/42-
Influenza B	Positive	18	0+/18-	16+/2-	0+/18-
SARS-CoV-2	Positive	102	0+/102-	0+/102-	99+/3-
SARS-CoV-2 + Influenza A/B	Negative	571	4+/567-	4+/567-	2+/569-

The NADAL® COVID-19 Ag+Influenza A/B plus Test showed no cross-reactivity between the parameters with the clinical specimens.

NADAL® COVID-19 Ag+Influenza A/B plus Test (COVID-19 Ag)	RT-PCR			
	Positive	Negative	Total	
	Positive	99	2	101
	Negative	3	629	632
Total	102	631	733	

Diagnostic sensitivity: 97.06% (91.71% - 98.99%)*

Diagnostic specificity: 99.68% (98.85% - 99.91%)*

Overall agreement: 99.32% (98.41% - 99.71%)*

*95% confidence interval

NADAL® COVID-19 Ag+Influenza A/B plus Test (Influenza A)	RT-PCR			
	Positive	Negative	Total	
	Positive	38	4	42
	Negative	4	687	691
Total	42	691	733	

Diagnostic sensitivity: 90.48% (77.93% - 96.23%)*

Diagnostic specificity: 99.42% (98.52% - 99.77%)*

Overall agreement: 98.91% (97.86% - 99.45%)*

*95% confidence interval

NADAL® COVID-19 Ag+Influenza A/B plus Test (Influenza B)	RT-PCR		
	Positive	Negative	Total
	Positive	16	4
Negative	2	711	713
Total	18	715	733

Diagnostic sensitivity: 88.89% (67.20% - 96.90%)*

Diagnostic specificity: 99.44% (98.57% - 99.78%)*

Overall agreement: 99.18% (98.23% - 99.45%)*

*95% confidence interval

Analytical performance

Analytical sensitivity study (detection limit)

The 95% detection limit of the NADAL® COVID-19 Ag+Influenza A/B plus Test for SARS-CoV-2 is 75.5 TCID₅₀/mL and was determined via dilution of inactivated SARS-CoV-2 virus.

Concentrations determined as the 95% detection limit for each influenza A/B strain tested are listed below:

Influenza A: Brisbane/59/2007 (H1N1), 1.32 × 10³ TCID₅₀/mL

Influenza A: Hong Kong/4801/2014 (H3N2), 5.68 × 10² TCID₅₀/mL

Influenza B: Brisbane/60/2008, 1.54 × 10³ TCID₅₀/mL

Influenza B: Phuket/3073/2013, 2.18 × 10³ TCID₅₀/mL

Measuring range

No prozone effect was observed when testing inactivated SARS-CoV-2 viruses (up to 1.51 × 10⁶ TCID₅₀/mL) and the following inactivated Influenza strains: H1N1 (2.64 × 10⁶ TCID₅₀/mL), H3N2 (1.136 × 10⁶ TCID₅₀/mL), B/Brisbane/60/2008 (3.08 × 10⁶ TCID₅₀/mL), B/Phuket/3073/2013 (1.136 × 10⁶ TCID₅₀/mL).

Analytical reactivity

The human, avian or other animal-derived influenza A and B virus strains listed below have been tested and showed a positive reaction with the NADAL® COVID-19 Ag+Influenza A/B plus Test.

Strains
SARS-CoV-2
A/Brisbane/59/2007 (H1N1)
A/New Caledonia/20/1999 (H1N1)
A/PR/8/34 (H1N1)
A/Solomon Islands/03/2006 (H1N1)
A/California/07/2009 (H1N1pdm)
A/Canada/6294/2009 (H1N1pdm)
A/Brisbane/10/2007 (H3N2)
A/HK/8/68 (H3N2)
A/Texas/50/2012 (H3N2)
A/Victoria/361/2011 (H3N2)
A/Wisconsin/67/2005 (H3N2)
A/Hong Kong/4801/2014 (H3N2)
A/Cambodia/R0405050/2007(08/216) (H5N1)
A/Netherlands/12/2000(07/336)
A/Hongkong/1073/1999(08/208)
A/Guangdong-Maonan/SWL1536/2019
A/Hong Kong/2671/2019
B/Brisbane/60/2008

Strains
B/Malaysia/2506/2004
B/Massachusetts/2/2012
B/Phuket/3073/2013
B/Texas/6/2011
B/Wisconsin/1/2010
B/Yamagata/16/1988
B/Washington/02/2019

Interfering substances

The following substances, normally present in respiratory specimens or artificially introduced into the respiratory tract, were evaluated at the concentrations listed below. All substances were mixed with a 1/20000 SARS-CoV-2 virus dilution (1.51 × 10⁶ TCID₅₀/mL) and showed no interference with the NADAL® COVID-19 Ag+Influenza A/B plus Test.

Substance	Concentration
Zanamivir	5 mg/mL
Oseltamivir	10 mg/mL
Artemether/Lumefantrine	50 µM
Doxycycline hyclate	70 µM
Quinine	150 µM
Lamivudine	1 mg/mL
Ribavirin	1 mg/mL
Daclatasvir	1 mg/mL
Mucin from bovine submaxillary glands, type I-S	100 µg/mL
EDTA-anticoagulated human blood	5% (v/v)
Biotin	100 µg/mL
Neosynephrin® (phenylephrine)	10% (v/v)
Afrin® Nasal Spray (oxymetazoline)	10% (v/v)
Saline nasal spray	10% (v/v)
homeopathic Zicam® Allergy Relief Nasal Gel	5% (v/v)
Sodium cromoglycate	20 mg/mL
Olopatadine hydrochloride	10 mg/mL
Paracetamol	199 µM
Acetylsalicylic acid	3.62 mM
Ibuprofen	2.425 mM
Mupirocin	10 mg/mL
Tobramycin	5 µg/mL
Erythromycin	81.6 µM
Ciprofloxacin	30.2 µM

Cross-reactivity

Specimens spiked with the following pathogens were tested using the NADAL® COVID-19 Ag+Influenza A/B plus Test:

Pathogen	Concentration
Influenza A H1N1 virus	6×10 ⁵ TCID ₅₀ /mL
Influenza A H3N2 virus	9×10 ⁵ TCID ₅₀ /mL
Influenza A H5N1 virus	8×10 ⁴ TCID ₅₀ /mL
Influenza A H7N9 virus	9×10 ⁴ TCID ₅₀ /mL
Influenza B	3×10 ⁵ TCID ₅₀ /mL
Adenovirus type 1	3×10 ⁵ TCID ₅₀ /mL
Adenovirus type 2	2×10 ⁵ TCID ₅₀ /mL
Adenovirus type 3	1×10 ⁶ TCID ₅₀ /mL
Adenovirus type 5	4×10 ⁵ TCID ₅₀ /mL
Adenovirus type 7	1×10 ⁶ TCID ₅₀ /mL
Adenovirus type 55	4×10 ⁵ TCID ₅₀ /mL
Respiratory syncytial virus type A/B	3×10 ⁵ TCID ₅₀ /mL
Human coronavirus 229E	3×10 ⁴ TCID ₅₀ /mL
Human coronavirus OC43	1×10 ⁵ TCID ₅₀ /mL

Pathogen	Concentration
Human coronavirus NL63	1×10 ⁴ TCID ₅₀ /mL
Human coronavirus HKU1	1×10 ³ ng/mL
MERS-Coronavirus Florida/USA-2/Saudi Arabia.2014	4×10 ⁴ TCID ₅₀ /mL
Parainfluenza virus type 1/2/3/4	1×10 ⁵ TCID ₅₀ /mL
Rhinovirus type A16	1×10 ⁵ TCID ₅₀ /mL
<i>Legionella pneumophila</i> Bloomington-2	1×10 ⁵ cells/mL
<i>Legionella pneumophila</i> 82A3105	1×10 ⁵ cells/mL
<i>Mycobacterium tuberculosis</i> K	1×10 ⁵ cells/mL
<i>Mycobacterium tuberculosis</i> Erdman	1×10 ⁵ cells/mL
<i>Mycobacterium tuberculosis</i> HN878	1×10 ⁵ cells/mL
<i>Mycobacterium tuberculosis</i> CDC1551	1×10 ⁵ cells/mL
<i>Mycobacterium tuberculosis</i> H37Rv	1×10 ⁵ cells/mL
<i>Streptococcus pneumoniae</i> 475298 [Maryland (D1) 6B-17]	1×10 ⁵ cells/mL
<i>Streptococcus pneumoniae</i> 178 [Poland 23F-16]	1×10 ⁵ cells/mL
<i>Streptococcus pneumoniae</i> 262 [CIP 104340]	1×10 ⁵ cells/mL
<i>Streptococcus pneumoniae</i> Slovakia 14-10 [29055]	1×10 ⁵ cells/mL
<i>Streptococcus pyogenes</i> T1	1×10 ⁵ cells/mL
<i>Mycoplasma pneumoniae</i> mutant 22	1×10 ⁵ cells/mL
<i>Mycoplasma pneumoniae</i> FH strain of Eaton Agent	1×10 ⁵ cells/mL
<i>Mycoplasma pneumoniae</i> M129-B7	1×10 ⁵ cells/mL

Precision

Repeatability

Repeatability was established by testing 10 replicates of negative, low and medium positive specimens using NADAL® COVID-19 Ag+Influenza A/B plus test lot.

Reproducibility

Reproducibility was established by testing triplicates of negative and high negative, as well as low, medium and high positive specimens. Testing was performed by 6 operators using 3 independent NADAL® COVID-19 Ag+Influenza A/B plus test lots at 3 different sites on 5 separate days.

The NADAL® COVID-19 Ag+Influenza A/B plus Test demonstrated acceptable repeatability and reproducibility. The negative and positive values were correctly identified >99% of the time.

15. References

- Cui J, Li F, Shi ZL, Origin and evolution of pathogenic coronaviruses, *Nat Rev Microbiol* 2019; 17:181-192.
- Su S, Wong G, Shi W, et al, Epidemiology, genetic recombination, and pathogenesis of coronaviruses, *Trends Microbiol* 2016;24:490-502.
- Weiss SR, Leibowitz JL, Coronavirus pathogenesis, *Adv Virus Res* 2011; 81:85-164.
- Murphy, B.R., and R.G. Webster, 1996, Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
- Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
- CDC website: <http://www.cdc.gov/flu/>
- Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. *The New England Journal of Medicine*, 353 (13):1363-1373.

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a) 20 hisopos estériles, CE 0197

Jiangsu Changfeng Medical Industry Co., Ltd
Touqiao Town, Guangling District, Yangzhou,
Jiangsu 225109 China (representante EU
autorizado Llins Service & Consulting GmbH, Obere
Seegasse 34/2, 69124 Heidelberg, Germany)

b) 20 hisopos estériles, CE 0197

Jiangsu Rongye Technology Co., LTD, Touqiao
Town, Yangzhou City, Jiangsu Province, China
(representante EU autorizado Riomavix S.L., Calle
de Almansa 55, 1D, Madrid 28039 Spain)

c) 20 hisopos estériles, CE 0197

CITOTEST LABWARE MANUFACTURING CO., LTD
No.48, Xinxu Road, Haimen, Jiangsu province
(representante EU autorizado WellKang Ltd,
Enterprise Hub, NW Business Complex,
1 Beraghmore Rd., Derry, BT48 8SE, Northern
Ireland)

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- a) 20 tamponi sterili, CE 0197

Jiangsu Changfeng Medical Industry Co., Ltd
Touqiao Town, Guangling District, Yangzhou,
Jiangsu 225109 China (rappresentante autorizzato
EU Llins Service & Consulting GmbH, Obere
Seegasse 34/2, 69124 Heidelberg, Germany)

- b) 20 tamponi sterili, CE 0197

Jiangsu Rongye Technology Co., LTD, Touqiao
Town, Yangzhou City, Jiangsu Province, China

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c) 20 tamponi sterili, CE 0197

CITOTEST LABWARE MANUFACTURING CO., LTD
No.48, Xinxiu Road, Haimen, Jiangsu province
(rappresentante autorizzato EU WellKang Ltd,
Enterprise Hub, NW Business Complex,
1 Beraghmore Rd., Derry, BT48 8SE, Northern
Ireland)

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15. Bibliografia

- Cui J, Li F, Shi ZL, Origin and evolution of pathogenic coronaviruses, Nat Rev Microbiol 2019; 17:181-192.
- Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses, Trends Microbiol 2016;24:490-502.
- Weiss SR, Leibowitz JL. Coronavirus pathogenesis, Adv Virus Res 2011; 81:85-164.
- Murphy, B.R., and R.G. Webster, 1996. Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
- Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
- CDC website: <http://www.cdc.gov/flu/>
- Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. The New England Journal of Medicine, 353 (13):1363-1373.

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15. Bibliografia

- Cui J, Li F, Shi ZL, Origin and evolution of pathogenic coronaviruses, Nat Rev Microbiol 2019; 17:181-192.
- Su S, Wong G, Shi W, et al, Epidemiology, genetic recombination, and pathogenesis of coronaviruses, Trends Microbiol 2016;24:490-502.
- Weiss SR, Leibowitz JL, Coronavirus pathogenesis, Adv Virus Res 2011; 81:85-164.
- Murphy, B.R., and R.G. Webster, 1996, Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
- Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
- CDC website: <http://www.cdc.gov/flu/>
- Anne Moscona. Neuraminidase inhibitors for Influenza, 2005. The New England Journal of Medicine, 353 (13):1363-1373.

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a) 20 zaragatoas estêreis, CE 0197

 Jiangsu Changfeng Medical Industry Co., Ltd
Touqiao Town, Guangling District, Yangzhou,
Jiangsu 225109 China (representante EU
autorizado Llins Service & Consulting GmbH, Obere
Seegasse 34/2, 69124 Heidelberg, Germany)

b) 20 zaragatoas estéreis, CE 0197

Jiangsu Rongye Technology Co., LTD, Touqiao Town, Yangzhou City, Jiangsu Province, China (representante EU autorizado Riomavix S.L., Calle de Almansa 55, 1D, Madrid 28039 Spain)

c) 20 zaragatoas estéreis, CE 0197

CITOTEST LABWARE MANUFACTURING CO., LTD No.48, Xinxiu Road, Haimen, Jiangsu province (representante EU autorizado WellKang Ltd, Enterprise Hub, NW Business Complex, 1 Beraghmore Rd., Derry, BT48 8SE, Northern Ireland)

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15. Lähteet

1. Cui J, Li F, Shi ZL, Origin and evolution of pathogenic coronaviruses, *Nat Rev Microbiol* 2019; 17:181-192.
2. Su S, Wong G, Shi W, et al, Epidemiology, genetic recombination, and pathogenesis of coronaviruses, *Trends Microbiol* 2016;24:490-502.
3. Weiss SR, Leibowitz JL, Coronavirus pathogenesis, *Adv Virus Res* 2011; 81:85-164.
4. Murphy, B.R., and R.G. Webster, 1996, Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
5. Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
6. CDC website: <http://www.cdc.gov/flu/>
7. Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. *The New England Journal of Medicine*, 353 (13):1363-1373.

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